
PicoPhi RCA Premix with Specific Primers

Cat#: PM100S or PM1000S



FOR RESEARCH USE ONLY

PicoPhi RCA Pre-mix Kit is intended for molecular biology use and *in vitro* use only. This product is not intended for diagnosis, prevention or treatment of a disease in human beings or animals.

Store Kit at <-20°C on Receipt

Description

This kit was developed to amplify specific template DNA using customized primers through Rolling Circle Amplification (RCA) technology. This premix is optimized for efficient and reproducible RCA, an equivalent of *PicoPhi* RCA with Specific Primers Kit (Cat#: Phi100S). This premixed formulation saves time and reduces contamination due to a reduced number of pipetting steps required for RCA set up. The starting circle DNA template concentration should be more than one nano-gram. Using excessive multiple primers complementary to template DNA in RCA reaction will generate multiple amplification origins and significantly enhance DNA amplification efficiency.

Kit Components:

Components	100 Reactions (Cat#: PM100)	1000 Reaction (Cat#: PM1000)	Storage
Sample Buffer	800 μ l	10x 800 μ l	-20°C
PicoPhi Premix with Specific Primers (#PM100S)	500 μ l	10x 500 μ l	<-20°C

Detailed Protocol

1. Preparation of Sample Mix:

Sample Mix could be prepared, depending on material sources, as described below:

1.1. Purified DNA or DNA ligation/assembly reactions:

Transfer 3 μ l of Sample Buffer into a 0.2 ml PCR tube.

Add 1 μ l of circular DNA (≥ 1 pg/ μ l) to the above PCR tube.

Heat to 95 °C for 3 minutes and then quickly cool to 4°C.

Keep the samples on ice until use.

1.2. Bacterial colonies:

Transfer 4 μ l of Sample Buffer into a 0.2 ml PCR tube.

Pick 1/10 to 1/100 of the colony (approximately $10^2 \sim 10^4$ cells) and add to the above PCR Tub.

Heat to 95 °C for 3 minutes and then quickly cool to 4°C.

Keep the samples on ice until use.

1.3. Liquid bacterial culture:

Transfer 3.0~3.8 μ l of Sample Buffer into a 0.2 ml PCR tube.

Add 0.2~1.0 μ l of saturated overnight culture to the above PCR tube.

Heat to 95 °C for 3 minutes and then quickly cool to 4°C.

Keep the samples on ice until use.

1.4. Glycerol stock:

Transfer 3.0~3.8 μ l of Sample Buffer into a 0.2 ml PCR tube.

Add 0.2~1.0 μ l of glycerol stock to the above PCR tube.

Heat to 95 °C for 3 minutes and then quickly cool to 4°C.

Keep the samples on ice until use.

Note: Heating at higher temperature or longer time may increase the probability of nicking target DNA and releasing host genomic DNA into cell lysis to the reaction, where the host genomic DNA will compete with the desired template DNA during amplification.

2. DNA Amplification:

Add 5 μ l of PicoPhi Premix (#PM100S) and 1 μ l of specific primer(s) to the above PCR tube contains 4 μ l of Sample Mix as showed below:

Component	Volume/reaction
Sample Mix with denature DNA	4 μ l
PicoPhi Premix with Specific Primers (#PM100S)	5 μ l
Specific Primer (100 μ M)	1 μ l
Final Volume	10 μl

Mix well and incubate at 30 °C for 3~24hrs.

3. Inactivate PicoPhi: Incubating at 65 °C for 10 min, and then cool to 4 °C.

4. Perform Downstream Application:

- The amplified DNA can be directly used for the cycle sequencing reaction without purification;
- The amplified DNA can be directly used for DNA restriction enzyme digestion;
- An aliquot of the amplified DNA can be examined by agarose gel;
- Specially, the amplified RCA products as concatemers, maybe after diluted, can be used to directly transform *Bacillus*. Compared to the traditional *bacillus* transformation methods, the RCA1.0 products would give the highest transformation rate, which is of critical importance to *bacillus* gene cloning, expression and library construction.

5. FAQ and Troubleshooting: Please contact us at info@evomicsscience.com.

Limited Use License and Warranty

Use of this Product is subject to the following terms and conditions. Evomic Science LLC ("We") warrant that the Product meets the specifications described in this manual. If the terms and conditions are not acceptable, please return all components of Products with original receipt by the original buyer to Evomic Science LLC within 10 business days. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. No right or license to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This Product should be used in accordance with the NIH guidelines developed for recombinant DNA and genetic research. Evomic Science LLC disclaims any and all responsibility for injury or damage which may be caused by the failure of the buyer or any other person to use the Product in accordance with the terms and conditions outlined herein. If you should have any questions or concerns about any of our products, please contact us at info@evomicsscience.com.